

See “Magnified single-balloon enteroscopy in the diagnosis of intestinal follicular lymphoma: a case series” on page ____.

Supplementary Table 1. Specifications of SIF-Y0007

	SIF-Y0007	SIF-Q260
Optical system		
Field of view	Wide: 140°; tele 60°	140°
Direction of view	Forward viewing	Forward viewing
Depth of field	Wide: 7–100 mm, tele: 2–3.5 mm	5–100 mm
Magnification	×80	×20
Distal end		
Outer diameter	9.9 mm	9.2 mm
Insertion tube		
Outer diameter	9.2 mm	9.2 mm
High force transmission	Yes	No
Bending section		
Angulation range (U/D/R/L)	180°/180°/160°/160°	180°/180°/160°/160°
Passive bending	Yes	No
Working length	2,000 mm	2,000 mm
Total length	2,345 mm	2,345 mm
Instrumental channel		
Inner diameter	2.8 mm	2.8 mm
Overtube	ST-SB1	ST-SB1

U, up; D, down; R, right; L, left.